From: Caloggero, Dina (DPH)

Sent: Thursday, May 20, 2010 12:48 PM

To: Nassif, Julianne (DPH); Smole, Sandra (DPH); Clemmer, Jill (DPH); Kane01, Peter (DPH);

Jenner, Jennifer (DPH); Servizio, Paul (DPH); Jacobsen, Patricia (DPH); Rubin, Alan (DPH); Alfonso, Felipe (DPH); Salemi, Charles (DPH); Stiles, Tracy (DPH); Belanger, Peter (DPH); Holland, James (DPH); Borne, Alan (DPH); Gauthier, Cheryl (DPH); Chen, Karen (DPH); Konomi, Raimond (DPH); Kazianis, Arthur (DPH); Guillet, Jasmine (DPH); Elvin, Paul (DPH); Hennigan, Scott (DPH); Krumholz, Glenn (DPH); Wang, Xingtai

(DPH)

Cc:DiNatale, Margaret (DPH); Han, Linda (DPH)Subject:Thermometer Calibration Verification QA SOP 016

The new Thermometer Calibration Verification QA SOP 016 is now available on \\dph-sli-fs2\shared\Public\Quality \\Assurance\QA SOPs \\Assurance\QA

- Record each of the 3 readings for both the NIST and Test thermometers on the Thermometer Calibration Verification Worksheet (Attachment A).
- The thermometer calibration is considered to be verified if at least 2 of 3 readings of the TEST thermometer are within ±2°C of the corresponding NIST reading. Laboratories may apply more stringent tolerance criteria (for example, ±1°C) depending on the temperature accuracy required for a particular test or storage unit.

The worksheet was also modified to reflect these changes and is included in the SOP. Laboratory supervisors must ensure that staff aware of this change to the procedure. Training can be documented using the Training Documentation worksheet included in QA SOP 013 as Attachment #3.

Please call me if you have any questions. I appreciate everyone's patience in getting this SOP in place.

Thanks...Dina